

Cell Counting and Particle Analysis from Image Captured by using Microscope

Swapnali Patil¹, Mayuri Pawaskar², Varsha Sawant³, Kajal Godage⁴, Prof.B.A.Jadhawar⁵

B.E Student, Department of Computer Science & Engineering, DACOE, Karad^{1,2,3,4}

Professor, Department of Computer Science & Engineering, DACOE, Karad⁵

Abstract: Blood cell count plays very important role in medical diagnosis. Increment or decrement in the count of blood cell causes many diseases to occur in the human body. There are different techniques of blood cell counting which involves conventional as well as automatic techniques. The conventional method of manual counting under microscope is time consuming and yields inaccurate results[5]. A solution to this problem, to provide software-based cost effective and an efficient alternative in recognizing and analyzing blood cells, This paper presents the preliminary study of automatic blood cell counting based on digital image processing. The number of blood cell count that is RBC & WBC count is then may be use to diagnose the patient as well as detection of abnormalities. For this purpose, few pre- processing and post-processing techniques have been implemented on blood cells image in order to provide a much clearer and cleaner image.

Keywords: Automatic detection, Biomedical image processing, Blood cell analysis, Segmentation.

I.INTRODUCTION

The observation of blood samples by expert operators is -The biggest detraction of microscopy, namely its one of the diagnostic procedures available for the recognition of different diseases. The human visual inspection is tedious, lengthy and repetitive, and suffers from the presence of a non-standard precision: that is because it depends on the operator skills. These reasons have limited its statistical reliability. On the other side, the automated analysis by computer requires only one image and not a blood sample; for this, it turns out to be less expensive but at the same time more scrupulous in providing precise standards [1].

This project describes a preliminary study of developing a detection of blood cell count using microscopic blood sample images. Analyzing through images is very important as from images, diseases can be detected and diagnosed at earlier stage. From there, further actions like controlling, monitoring and prevention of diseases can be done.

II.NEED OF PROJECT

To detect and find a blood cell count and produce an accurate cell count report. This would be very helpful to a physician in identifying the cause of his patient's diseases. Analysis and processing of a microscopic image, in order to provide an automated procedure to support the medical activity and easily identify the intensity of diseases. Also maintain the patient detail and manage detail.

III. GOALS AND OBJECTIVES

A.Goals

-The main goal of this work is the analysis and processing of a microscopic image, in order to provide an automated procedure to support the medical activity and easily identify the intensity of diseases.

- dependence on the skill, experience and motivation of a human technician is to be removed. Used with an automated digital microscope, which would allow entire slides to be examined, it would allow the system to make diagnoses with a high degree of certainty.
- B. Objective
- -To provide this system is user friendly so that the pathologist can observe the blood cell sample image in different view. Also pathologist can get the result of blood cell test within seconds as tested.
- -To count the number of cell automatically and provide the accurate number of cells.

IV.RELATED WORK

The conventional method of manual counting under microscope is time consuming and yields inaccurate results. Manual method for counting WBC and RBC.

A. White Blood Cell

Calculation of the WBC count is achieved by following the formula below using these factors:

- The number of WBCs counted in the big squares i.
- The dilution of the cell solution ii.
- The number of counted big squares iv. The volume iii. above a big square

iv. The volume above a big square

WBC"s/µl=

Number of counting WBC's * dilution

Number of counted square * volume above one small square

B. Red Blood Cell

The calculation is achieved by following the formula below using these factors:



- i. The number of RBCs counted in the small squares
- ii. The dilution of the cell solution

iii. The number of counted small squares

iv. The volume above one small square

RBC"s/µl=

Number of counting BC's * Dilution

Number of counted square * volume above one small square

V. PROPOSED METHODOLOGY

A. Functionality of proposed methodology

The overall Methodology of image processing techniques in the proposed method is described in below.

First of all an input image of already prepared slide is take from the camera attached to the microscope, then through different image processing techniques and image filters the extra unwanted information is removed from that input image. After removal of noise the area of interest that is the area where we need to count the blood cells are cropped from the image and by labeling algorithm of connected pixels blood cells can be count to put in different equations to abstained cell count report.



1. Image Pre-processing :

Step 1: Input image

In this, when you click on browse button- the current directory is opened, and from this you can select any patient's folder to calculate number of blood cells in his body. Like each patient's folder consists of blood images with the medical report of that patient.

Step 2: RGB to Gray conversion

Image read in step-1 is in RGB format. This image is to be converted in to grayscale format. After converting this image in to grayscale, practically source of noise is thermal noise and artifacts i.e. noise added due to body movement.

Step 3: Median Filtering

Median filtering is used to remove noise which is added at the time of capturing the microscopic blood images. This unwanted noise can be the dust particles that are present on the blood slide while preparing the slides.

2. Image Segmentation :

Step 1: Gray Thresholding

In this using thresholding algorithm removes the unwanted cells and extract and segment the image object and also clean the image [5].

Step 2: Filling Holes

Holes presented in the binary image form. This image is filled so that the blood cells are segmented properly.

3. Image Post-processing:

Step 1: Removing Borders

The cells occupying the border contains less information and has to be removed to reduce the complexity. The segmented image is labeled.

Step 2: Image Labeling

Labels the connected objects. The objects can have a value of either 4 or8, where 4 specifies 4-connected objects and 8 specifies 8-connectedobjects.

B. Modules and Their Functionalities:

Modules:

This system consists of three modules those are:

- A. Patient Information Management 1. Load Patient Information
 - 2. Set Test of patient.
- B .Update information of Patient 1 .Blood Cell Counting
- C. Labeling and Particle Analysis. 1. Printing

A. Patient Information Management



Fig. 2. Login page of the System



1. Patient Details:

5 5 File Patient1	Master Partical Size Analysis	Microlmage Processing Settings Help	Logout	- a ×
Options New Patien Patient Name	et O Search Patient	Tet Deals		Ad fee
Patient Name Addees Contact No. Gender Releved By Basination Rec. Specimen Tachinojan Date Date Patiet Te	Mayor A P-Ne 16 Face De Sare 21 1900/1868 Face v Inc v		E Teo piter cust	
Copyright © SAG	R.O. Soft (SR-11 Model, Indian D	esign Patent No 263018)		Calibrated at 100X - 3/100/100 pixel - 1/100/100 pixel
- 6	- 🔛 📋 🔹	o (5)		• 2 10 (C 00.6

Fig.3 Patient Details

In this module add information about patients for e.g. add patient name, address, and age, contact no, gender, referred by examination req, specimen, technician, and pathologist. After that save these information.

Suppose we want to search any particular patient information then you can easily search the patient. In this module you can add test, in this display the test details like group name, test description resultant count, normal range and comment.

In group name various test are included like W.B.C. count report, Haematology, Urine Examination Report. In add test we have get image.

In get image we can browse image select the camera then capture the image. Here display the cell details, cell range details, print cell details, print image, process in MS PAINT, Export image, Finalize result.

2. Print Patient List

In these generate the report and print the patient list.

B.Test Management: Test Groups

The Mentality Reschicable Machine Mentality Mentality of the Lapid The Section Terminal Section of the Control Terminal Secti			- 0
Instruction Set Mark Text Sectors Set Mark Text Sectors NOW Implementation Set Mark Text Sectors Set Desce Implementation Set Sectors Set Desce Implementation Set Sectors Set Desce Set Desce Set Sectors			
Variat Say 1	ie Patient Master Partical Size Analysis	Micro Image Processing Settings Help Logout	
view w w w w w w w w w w w w w w w w w w	- 144		
a Deser b Deser	e onde	H Grandberteten	
Image: Section of the Control of th	a Description Replaced	DIFFERENTIAL COUNT	
image:	transfer. Land	2. SMEAR EXAMINATION	
india Construction With india Construction With india Construction india Construction indindia Construction indindia Con		CINY YOOJOTAMINH C	
Bre Open Profile Average Profile Average		4 URINE EXAMINATION REP.	
Bite Open PUT MARAY 10 PUT MARAY 10 Activate Windows Solid bit 50-11Mee wide Texplorents		5 STOOL ANALYSIS REPORT	
See One protocology Activate Windows Cell of Antipal Activate Windows Cell of Antipal Activate Windows Cell of Antipal Activate PP Statust Mindows Activate Windows		6 WBCCOUNT REPORT	
Activate Windows Gala X comparative get 5 50x0 for 50x 11 Mere Non Segnitizen No. 30x10 Celleret # 110x - 1000 (100 - 100 Mere)	Save Close	7 SPUTAMANALYSIS	
Activate Window. Ga IX, https://www.google.com/google.c			
AcTuats Windows Data 75,500,516,500,1104ed New Supp New Yor - 2010 Collect at the Today State St			
Activate Window. Gai X, Endoge Section 2010 get \$ \$60.01x150x1144ex New SuppherNs .3010 Califord TB: 10.1010get .10010g			
Activate Windows Socials & Company & Activate Socials & Company & Activate Socials & Company & Statute Colleger # 105-1106/jpg-1000100			
Activate Windows Grup 17 writige a skint get 2 MAD lef 50, 11 More New Despitaertis - 20110 Galande et 100 - 1400 Higher - 1400 Higher			
Activate Windows Color & Control and Activate Windows Set 5 500 Set 50. Thirde New Registrate No. 2010 Colored TEXT-NUMBER - 1990 TR			
Activate Westware Gen PC writing in a bein giel SAD bef 50, 1166et blan being Raertin - 50010 Galande et 100 - 50010 part - 100100 part			
Activate Windows Only 75 Status (Section 2001) Galaxies (Section 2001) Galaxies (Section 2001)			
Activate Windows On VP writing in active pet 9.400.046 500,11164ek biom.beografikeet No30010 Galiwel et Wo-300100 pair -300100			
Activate Windows Gala X company and the State St			
grit © 540,0 Set (50,-11 Mode), Indian Design Patent No 200410 prest - 9130/100 prest - 9130/100 pr			
	ight © SAGLO Soft (SGL-11 Model, Indian)	Design Patent No 200510	Calibrated et 100X - X:100/100 pixel - V:100/100 pixel

Fig .4 Group Detail

In these we can add group description like HAEMATOLOGY, W.B.C. Count Report, SPUTAM Analysis and save the group description.

3. Test Information

In test information add group name, Test description Result with unit, Normal Range and save the test information.

Options O New Patie	t 🖲 Seach Patient	Test Details									
Patient Name 10	nut v	100,0	Group_hane	Test, Name	Cogonal (mage	Processed	Joage read	UN	Converts	Lifey_Late	
	Seath										
Dataset Information					Assign Te	1st		1			
Patient Name	Mayuri			Test Details							
Address	A/P-Vite, Tai-Patan, Dati-Satara			Group Name	DIFFERENTIAL CON	UNT		1			
				Tell Description				1			
Age Contract No.	21			Resultant Count	a unit		Gatimanal				
Center	Example V			Normal Parcen			Count				
Edward By	nfastantenda v										
Comparison Days	81			Connerts							
Spectreen	87										
Technician	Im				New Test	Save Te	st Close				
Pathologial	(nd			_	_		_				
Date	16-10-2015 04:46:02										
Dist											
PHETE	z bare Clove										
			_	_	_	_	_	_			_

Fig.5 Add Test

Organization Details

				- 0
File Patient Master Partical Size Analysis	Micro Image Processing Settin	ings Help Logout		
		Oronization In	formation	
		Ormitation Data		-
		Name	DALDE	
		Address	kand	
		Ovner Nane	abcd	
			Save Close	
opyright © SAGLO Soft (SGL-11 Model, Indian D	esign Patent No 263618)			Calibrated at: 100X - X:100/100 pixel - Y:100/100 pixel
- 🤗 🕨 🗀 🛛	a 🕾 .			

Fig.6 Organization Details

In Organization details add name of Organization, Address, Owner name and save the organization details.

Application User

Application User					0	
Der Honaton Full Name	Neur Permiter	Useria	и И 3	held	Nave	
User Name	Marie		2. P	KGLO	SAGLO Recently	
Passort	-					
THE DIST FRANCE	Save Close					

Fig. 7 Application User

In Application user add User information like a Full name Password, Re-Enter Password and save the application user.

3.Cell Counting

2	Image Processir
Browse Image	
Camera Settings Select Camera	
Integrated Webcam	
Start Capture	South State
Process Image	
Minimum Radius (1)	Sola Sac a a a gala a gala
	dealers and the second s
Maximum Radus (15)	
Result	
Cell Details	
Cell Range Details	
Print Cell Details	
Print Image	
Process in MS Paint	
Export Image	
Finalize Result	

Fig. 8 Capture Image



In this capture or browse the image and count the cell

VI.APPLICATIONS

This software used in Pathology Lab, Biomedical Center, Pharmacy College, Hospitals to count the number of cell automatically and provide the accurate number of cells and also manage the patient details management.

VII.FUTURE WORK

In this work it has been proposed an innovative method for a completely automatic identification of leukocytes from microscopic images, in order to provide an automated procedure as support for medical activity. The results obtained show that the proposed method is able to identify in a robust way the blood cells present in the image.

Proposed method for cell counting blood cells is similar for counting remaining two types of cell counting and through the same procedure we can develop another algorithm for detecting malaria parasite in the blood, malaria parasite also found in blood and the current method of detecting malaria parasite is to prepare the slide and study is under microscope to recognize malaria parasite in the blood, by improving our proposed algorithm we can also detect malaria parasite through computer vision technologies and techniques of image processing techniques. But this procedure would need a little advance microscope then the ordinary used in the clinical labs to achieve 100% accurate results.

VIII.CONCLUSION

This system has implemented the segmentation part and after testing this system it is concluded that, it is more timely efficient than the existing systems. It is reliable and cost effective than automated method. The aim of this system is to provide CBC i.e. complete blood Cell Count, which has been achieved by using the Powerful Image processing technique.

This system is more efficient in reducing the valuable time than manual system. The system is user friendly so that the pathologist can observe the blood cell sample image in different views. Also pathologist can get the result of blood cell test within 20 seconds to 1 min. as tested.

REFERENCES

- LorenzoPutzu, and Cecilia Di Ruberto, "White Blood Cells Identification and Counting from Microscopic Blood Image", 2013.
- [2] Sanaullah Khan, Faisal Saleh Khattak, Arslan Naseem, "An Accurate and Cost Effective Approach to Blood Cell Count", July 2012
- [3] Pallavi T. Suradkar ,"Detection of Malarial Parasite in Blood Using Image Processing ",April 2013.
- [4] Pooja R. Patil1, G.S. Sable2, Gauri Anandgaonkar3,"Counting of WBCs and RBCs from Blood Images using GrayThresholding", April 2014.
- [5] Sumeet Chourasiya,G.Usharani "Automatic /red Blood Cell Counting using water shead segmentation"2014.www.google.com